




Elecsys[®] Phospho-Tau (217P) Plasma

The power of ONE assay.
Redefining Alzheimer's diagnosis.



Elecsys® Phospho-Tau (217P) Plasma

The power of ONE assay: standalone, real-world clinically validated, IVD-cleared blood test with a unified dual-cut-off designed to rule-in and rule-out Alzheimer's-associated amyloid pathology across settings.^{1,2}

Alzheimer's disease (AD) is a progressive, neurodegenerative disease associated with cognitive dysfunction and behavioural impairments. It represents the most common form of dementia, contributing to at least 60 % of cases.² The accumulation of β -amyloid and tau are regarded as one of the earliest signs of the AD pathological cascade, taking place decades earlier prior to symptom onset.³⁻⁵

It has been demonstrated that Tau phosphorylation at threonine 217 occurs as a reaction to amyloid- β plaques and thus is closely linked to amyloid- β pathology.⁶ Specifically, it has been shown that pTau217 levels rise closer to the time when amyloid- β plaques become detectable by positron emission tomography (PET) imaging than to the time when neurofibrillary tangles (NFT) are detectable.⁶⁻⁸ Therefore, pTau217 has been proposed as a marker of amyloid- β proteinopathy.⁹

The pTau217 blood-based biomarker is considered a core one diagnostic biomarker in the Alzheimer's Association research framework¹⁰; the International Working Group considers the pTau217 biomarker as one of the most specific for AD pathology.⁷

pTau217 is an accurate blood-based biomarker with strong potential to improve routine clinical diagnosis of AD, particularly as an initial screening or triage test before confirmatory CSF (Cerebrospinal Fluid) or PET testing.¹⁰

Elecsys® Phospho-Tau (217P) Plasma redefines the current diagnosis of Alzheimer's disease

Elecsys® pTau217 is the best performing plasma assay in terms of detection of amyloid pathology (ROC AUC analysis)¹¹



Accurate and robust

- Clinically validated, IVD-cleared, real-world / routine samples
- High performance across diverse populations



Simplicity across settings

- Simple "draw and go" sample handling without pre-analytical limitations
- Standalone assay with unified dual cut-off across primary and secondary care settings



Enhanced clarity


- Allowing targeted CSF and PET testing to effectively identify amyloid pathology⁷⁻⁹
- Supporting diagnosis of early stages of AD for improved treatment options

Elecsys® Phospho-Tau (217P) Plasma redefines the current diagnosis of Alzheimer's disease¹

Ensuring demographic equity and unrestricted access by providing an operationally simple test that works consistently across diverse populations and care-settings – removing diagnostic access barriers.


958 total participants

Amyloid-PET positive
(n = 378)
Amyloid-PET negative
(n = 580)




Gender

Male
(n = 431)
Female
(n = 527)




eGFR

≥ 60
(n = 824)
45–59
(n = 91)
30–44
(n = 27)
≤ 29
(n = 1)
Missing
(n = 15)




Race

White
(n = 811)
Black or African American
(n = 110)
Asian
(n = 15)
Other
(n = 21)
Missing
(n = 1)




Clinical stage assesment (based on MMSE)

MMSE 28–30
(n = 295)
MMSE 25–27
(n = 363)
MMSE 22–24
(n = 227)
MMSE 19–22
(n = 73)




Care setting

Secondary/tertiary
(n = 720)
Primary
(n = 238)



Ethnicity

Not Hispanic or Latino
(n = 682)
Hispanic or Latino
(n = 188)
Not reported
(n = 5)
Unknown
(n = 1)
Missing
(n = 82)



Body mass index (BMI) kg/m²

Mean (SD)	27.9 (5.55)
Median	27.0
Q1–Q3	24.2–30.6
Min–Max	16.0–61.4
Missing	n = 72 (7.52%)

A unique approach to validation of Elecsys® Phospho-Tau (217P) Plasma

Clinical performance validated in diverse populations that match characteristics of real patients



Broad population of symptomatic patients evaluated for cognitive decline.¹



Included patients that are normally excluded in research cohorts.¹



30 sites (18 of which were used for the validation of the data) collected the intended use population (race heterogeneity, multiple comorbidities).¹

	Primary Care (low prevalence of 21.4%)	Secondary Care (moderate prevalence of 45.4%)	
Negative predictive value (NPV)	98.6%	93.7%	← Prevalence of amyloid pathology (pretest probability) affects NPV and PPV performance
Positive predictive value (PPV)	78.3%	90.2%	←
Positive percentage agreement (PPA, "sensitivity")	96.1%	94.2%	
Negative percentage agreement (NPA, "specificity")	94.7%	93.1%	


Primary Care patients

Subject self-refers/presents to the site or was referred by a primary care physician and has not seen a specialist before.

Secondary/Tertiary Care patients

Subject presents to centre or is recruited from a patient repository or has been referred by a specialist physician or has seen a specialist (neurologist, psychiatrist, geriatrician) before.

Based on the recently published CEOi Performance Standards for Blood-based Biomarkers (BBM, below), Elecsys® pTau217 is fulfilling all proposed criteria for triaging patients across care settings to confirmatory testing¹²





Accurate and robust blood test for patient management¹

AUC (%) of plasma biomarkers with increasing added Coefficient of Variation (CV).

All plasma biomarkers except for plasma Aβ42/Aβ40 showed less than 10% decline in AUC at 15% added CV and were thus considered stable.





Figure: presents AUC of plasma biomarkers with increasing added Coefficient of Variation (CV). All plasma biomarkers except for plasma Aβ42/Aβ40 showed less than 10% decline in AUC at 15% added CV and were thus considered stable.

Elecsys® pTau217 – A biomarker designed to support early diagnostic decisions for those who benefit most from therapeutic interventions¹³⁻¹⁵

Stage	Preclinical AD		MCI due to AD	Mild AD Dementia	Moderate AD Dementia	Severe AD Dementia
NIA-AA staging	Stage 1	Stage 2	Stage 3	Stage 4	Stage 5	Stage 6
MMSE score	30–27		27–25	24–19	18–10	9–0
No. of patients	295		363	300		
						

AUC: Area Under the Curve; GFAP: Glial Fibrillary Acidic Protein; NFL: Neurofilament Light Chain; MCI: Mild Cognitive Impairment; MMSE: Mini-Mental State Examination; NIA-AA: National Institute on Aging and Alzheimer’s Association

References

1. Elecsys® phospho-Tau (217P) plasma method sheet Roche Diagnostics GmbH.
2. Brum WS, et al. *Nat Aging*. 2023;3:1079-90.
3. Chong FP, et al. Tau Proteins and Tauopathies in Alzheimer's Disease. *Cell Mol Neurobiol*. 2018 Jul;38(5):965-980.
4. Smet-Nocca C. 1st Ed. *Tau Protein: Methods and Protocols; Methods in Molecular Biology*; Springer New York: New York, NY, 2017; Vol. 1523.
5. Wang Y, Mandelkow E. Tau in physiology and pathology. *Nat Rev Neurosci*. 2016 Jan;17(1):5-21.
6. Pichet Binette A, et al. Amyloid-associated increases in soluble tau relate to tau aggregation rates and cognitive decline in early Alzheimer's disease. *Nat Commun*. 2022 Nov 4;13(1):6635.
7. Barthélemy NR, et al. A soluble phosphorylated tau signature links tau, amyloid and the evolution of stages of dominantly inherited Alzheimer's disease. *Nat Med*. 2020 Mar;26(3):398-407.
8. Mattsson-Carlgren N, et al. Aβ deposition is associated with increases in soluble and phosphorylated tau that precede a positive Tau PET in Alzheimer's disease. *Sci Adv*. 2020 Apr 15;6(16):eaaz2387.
9. Theriault J, et al. Association of Phosphorylated Tau Biomarkers With Amyloid Positron Emission Tomography vs Tau Positron Emission Tomography. *JAMA Neurol*. 2023 Feb 1;80(2):188-199.
10. Jack CR Jr, et al. Revised criteria for diagnosis and staging of Alzheimer's disease: Alzheimer's Association Workgroup. *Alzheimers Dement*. 2024 Aug;20(8):5143-5169.
11. Ashton, N. J., Keshavan, A., Brum, W. S., et al. (2025). The Alzheimer's Association Global Biomarker Standardization Consortium (GBSC) plasma phospho-tau Round Robin study. *Alzheimer's & Dementia*. <https://doi.org/10.1002/alz.14508>
12. Schindler, S. E., Galasko, D., Pereira, A. C., et al. (2024). Acceptable performance of blood biomarker tests of amyloid pathology – recommendations from the Global CEO Initiative on Alzheimer's Disease. *Nature Reviews Neurology*, 20, 426–439. <https://doi.org/10.1038/s41582-024-00977-5>
13. Ding Z, Leung PY, Lee TL, Chan AS. Effectiveness of lifestyle medicine on cognitive functions in mild cognitive impairments and dementia: A systematic review on randomized controlled trials. *Ageing Res Rev*. 2023 Apr;86:101886. doi: 10.1016/j.arr.2023.101886. Epub 2023 Feb 18. PMID: 36806378
14. Barnett JH, Lewis L, Blackwell AD, Taylor M. Early intervention in Alzheimer's disease: a health economic study of the effects of diagnostic timing. *BMC Neurol*. 2014 May 7;14:101. doi: 10.1186/1471-2377-14-101. PMID: 24885474; PMCID: PMC4032565.
15. Rafii MS, Aisen PS. Detection and treatment of Alzheimer's disease in its preclinical stage. *Nat Aging*. 2023 May;3(5):520-531. doi: 10.1038/s43587-023-00410-4. Epub 2023 May 18. PMID: 37202518; PMCID: PMC11110912.

THESE PRODUCTS ARE NOT AVAILABLE FOR PURCHASE BY THE GENERAL PUBLIC.
ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS FOR USE.

COBAS and ELECSYS are trademarks
of Roche.

© 2026 Roche

MC-AU-04748 | June 2026

Published by:

Roche Diagnostics Australia Pty Limited
2 Julius Avenue North Ryde NSW 2113
Ph: +61 2 9860 2222
ABN 29 003 001 205

diagnostics.roche.com